

# DESIGN AND SYNTHESIS OF THE COMPREHENSIVE FRAGMENT LIBRARY A 3D ENABLED LIBRARY FOR MEDICINAL CHEMISTRY DISCOVERY

Warren S Wade<sup>1</sup>, Todd Meyer<sup>1</sup>, Peter Pallai<sup>1</sup>, Paolo Tosco<sup>2</sup>, James Zapf<sup>3</sup>, Laura Lingardo<sup>3</sup>, Gordon Alton<sup>3</sup>

<sup>1</sup>BioBlocks, Inc., <sup>2</sup>Cresset, <sup>3</sup>Visionary Pharmaceuticals

**BioBlocks**  
A CHEMISTRY PARTNER IN DRUG DISCOVERY

## INTRODUCTION



The Comprehensive Fragment Library (CFL) is a set of small, rigid, medicinally interesting fragments. This library originates from a starting set of >3 million potentially synthesizable virtual fragments designed from first principles and 3D enabled to maximize exploration of target interactions. Extensive clustering analysis allows broad coverage of chemistry space and provides an immediate follow-up strategy from any screening hit.

Supplementing Screening Decks with 3D compounds is an enormous task

Even common 2D cores are distributed unevenly in the literature

AbbVie Screening Collection, 2007<sup>1</sup>

193k Monosubstituted Quinolines in SciFinder



## BUILDING THE VIRTUAL FRAGMENT SET

An exhaustive set of ring structures containing ≤18 heavy atoms with 1 handle atom was built in BIOVIA Pipeline Pilot with several inclusion criteria:

- ≥1 ring, including bridged, spiro and fused ring connections
- ≤3 unique rings
- ≤2 rotatable bonds
- C,N,S,O atoms only
- ≤2 ring assemblies
- 1 atom has all rotatable bonds
- ≤1 S

**Rigid, partially aromatic structures are the highest value 3D enabling subset**

Core structures decorated with handles chosen for synthesis potential

| High Value          | Moderate Value   | Follow up only                 | Not chosen                |
|---------------------|------------------|--------------------------------|---------------------------|
|                     |                  |                                |                           |
|                     |                  |                                |                           |
| Neutral, extendable | Halogen, high MW | Interaction with handle likely | Modification loses charge |

The raw set of > 50 million structures with Me handles was filtered by medicinal chemistry criteria to produce a 7 million candidate fragment set

| Example CFL Fragments | Excluded | Reason                    |
|-----------------------|----------|---------------------------|
|                       |          | Two handles*              |
|                       |          | Unstable                  |
|                       |          | Isolated double bond*     |
|                       |          | Too many rotatable bonds* |

\*Potential follow up for a CFL fragment hit

Practical computational limit is <1000k diastereomers


Requires >600 CPU-years for full 3D similarity matrix

A 580k subset was selected for clustering:

Reduced mean of 16 heavy atoms to 14.6


Includes representatives of all 9000 ring types

Includes 100% of <12 atom fragments to 1.5% of 18 atom fragments



## CLUSTERING

830k diastereomers were generated from the 580k 2D fragments  
RDKit was used to enumerate diastereomers, with manual correction  
Cresset's XedeX generates s5 conformations for each




2D and 3D clustering methods were investigated on a 20k distributed subset  
4k clusters using a parallel C++ implementation of the k-medoids/CLARANS algorithm<sup>3</sup>  
Cluster tightness and separation assessed by the silhouette metric<sup>4</sup>

1.0: ideal; -1.0: misclustered

2D method: ECFP4

3D method: 0.75 fields, 0.25 shape




## 3D clustering method

- Alignment to methyl handle:
- Rotate structures around the handle
- Determine maximum similarity
- Similar compounds:
- Have common vectors
- Represent alternative sprouting choices

## 3D clustering produced better clusters

Comparison of two 6-member clusters

Silhouette: 0.509



## 3D clustering of the entire set was still computationally expensive

Initial clusters from a full similarity matrix of 150k compounds 14 CPU-years


Examined 10k, 25k and 50k clusters

10k gave a well-distributed 150k compound set; not scalable to the entire library

50k gave cluster sizes too small to be stable

Remaining compounds placed in 25k clusters of most similar medoid 20 CPU-years


PseudoSilhouette calculated because full matrix unavailable



## BUILDING THE PHYSICAL LIBRARY


### Compounds prioritized by HBond pattern

CFL fragments with 1 likely interaction motif preferred  
HBonds most common interactions with different targets  
Important for selectivity  
2D similarity measure  
Same patterns are distributed over clusters  
~50% of clusters have simple, valuable patterns



### HBond patterns were generated for all compounds

Pipeline Pilot Molecular Pharmacophore Fingerprints Component  
HBA, HBD paired when separated by 1-6 bonds  
35 core patterns including Named Functional Groups



### Available Chemical Directory™ (ACD) Search


Handles of each CFL Candidate generated and searched individually in DiscoveryGate™  
Grouped by HBond pattern, starting with least complex: >850,000 queries  
Hits filtered for real vendors with a real price: 11,184 results (1.6%), 2,618 unique cores  
Manually reviewed for medicinal chemistry value and purchase

### Available Compounds found in 1395 clusters (5.6% of total)

905 clusters have only 1 commercial member, 81% are 3D enabled

Others evenly distributed 2D and 3D by PMI1,2

On rod/disk axis = slope from 2-butynol or benzene near -1




CFL contains the high value commercial set in <1 plate

Most CFL compounds are not commercial

Can be synthesized from available starting materials

Library members are novel to targets, novel for IP




## PLATE 1

### Representative compounds picked from available clusters

Compounds >0.8 similarity to medoid  
Ring systems varied over columns  
HBond patterns grouped in rows  
Plate assembly underway  
Range of vendors  
Synthesis when handle options limited  
Purities confirmed  
200 mM DMSO solubility confirmed


|               | Plate 1 | CFL Design |
|---------------|---------|------------|
| # Fragments   | 55      | ~1000      |
| # Atoms       | 11.5    | 13         |
| AlogP         | 1.0     | 1.8        |
| HBD           | 0.91    | 0.6        |
| HBA           | 1.80    | 2.2        |
| Rot Bonds     | 0.91    | 0.5        |
| %3D           | 70%     | 80%        |
| %Unique Rings | 55%     | 75%        |
| % Commercial  | 80%     | 15%        |
| % Handles     | 100%    | 100%       |



## SCREENING AND FOLLOWUP

Library screened at 2 mM in a high concentration kinase inhibition assay  
Compared to a screen of the previous commercial library

### Fragment Assay Summary (51 to date)




Lower hit rate from added topology

Higher % confirm in dose response

Follow up possible on weak, interesting compounds


### Fragment Screen Summary (253 to date)



Privileged chemotypes present in both

### Representative Hit

Included as a positive control for kinases



### CFL clustering reveals many related compounds for rapid fragment SAR

| Independent Analogs  | CFL000025 | 10 Fragment Hits | Plate 1 |
|----------------------|-----------|------------------|---------|
| Handle Analogs       | 10        | 100              | 550     |
| Handle Vectors       | 4         | 47               | 257     |
| Same 2D Scaffold     | 37        | 343              | 3308    |
| Same HBond Pattern   | 274       | 12857            | 42667   |
| Same 3D Cluster      | 47        | 573              | 6371    |
| Total CFL Compounds  | 873       | 13676            | 50669   |
| % Unique to Compound | 100%      | 74%              | 26%     |
| CFL Analogs          | >30,000   | >1 M             | >25 M   |

### Leap-to-Lead™ generates fragment analogs

Compounds from the fragment hit column were assayed:

- New SAR identified from handle vectors
- New scaffolds identified from core alterations
- >3-fold enrichment of actives

~10-fold increase in best potency, increase in LE

Analogs include non-commercial compounds

IP identified despite well-known cores

### Fragment Assay Summary (89 to date)



One series was progressed to medicinal chemistry

Advanced to achieve μM potency in cell assays

Provisional filed, synthesis ongoing

## CONCLUSIONS

Libraries of commercial compounds are insufficient to generate new chemical matter  
The Comprehensive Fragment Library gives greater access to 3D-enabled interaction space  
Preliminary CFL screening generates high quality hits through its coverage of chemistry space  
Enhanced 3D enablement reveals new, unexpected chemotypes for important drug targets  
The CFL gives high quality entry points to our proprietary Leap-to-Lead™ discovery platform

## REFERENCES

- Akritopoulou-Zanze, I.; Metz, J.T.; Djuric, S.W. *Drug Discovery Today*, 2007, 948-952.
- Metz, J.T. Principal Moments of Inertia Protocol, Pipeline Pilot Community.
- Ng, R. T.; Han, J. *IEEE T. Knowl. Data. En.* 2002, 14, 1003-1016.
- Rousseuw, P. J. *J. Comput. Appl. Math.* 1987, 20, 53-65.

**CONTACT** **BioBlocks**  
A CHEMISTRY PARTNER IN DRUG DISCOVERY

BioBlocks is developing the CFL for use in Lead Discovery Collaborations. While the library is not available for independent purchase, we welcome new collaborators. For additional information or to discuss using the CFL and the Leap-to-Lead™ platform in a drug discovery effort, please contact [wwade@bioblocks.com](mailto:wwade@bioblocks.com) or visit our website: [www.bioblocks.com](http://www.bioblocks.com)